

## Product datasheet for **GA103370**

### Human LIS1 (PAFAH1B1) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PAFAH1B1 CRISPRa kit - CRISPR gene activation of human platelet activating factor acetylhydrolase 1b regulatory subunit 1
Symbol:	PAFAH1B1
UniProt ID:	<a href="#">P43034</a>
Synonyms:	LIS1; LIS2; MDCR; MDS; NudF; PAFAH
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)