

Product datasheet for **GA103375**

Human PAH activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PAH CRISPRa kit - CRISPR gene activation of human phenylalanine hydroxylase |
| Symbol: | PAH |
| UniProt ID: | <u>P00439</u> |
| Synonyms: | PH; PKU; PKU1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)