

Product datasheet for **GA103383**

Human Contactin 3 (CNTN3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CNTN3 CRISPRa kit - CRISPR gene activation of human contactin 3
Symbol:	CNTN3
UniProt ID:	Q9P232
Synonyms:	BIG-1; PANG; PCS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)