

## Product datasheet for **GA103384**

### Human **REG3A** activation kit by **CRISPRa**

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | REG3A CRISPRa kit - CRISPR gene activation of human regenerating family member 3 alpha |
| Symbol:       | REG3A  |
| UniProt ID:   | <a href="#">Q06141</a>   |
| Synonyms:     | HIP; HIP/PAP; INGAP; PAP; PAP-H; PAP1; PBCGF; REG-III; REG3                            |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                |

#### Product images:



[View online >](#)