

Product datasheet for **GA103386**

Human Parkin (PARK2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRKN CRISPRa kit - CRISPR gene activation of human parkin RBR E3 ubiquitin protein ligase
Symbol:	PRKN
UniProt ID:	O60260
Synonyms:	AR-JP; LPRS2; PDJ; PRKN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)