

## Product datasheet for **GA103391**

### Human PAX3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PAX3 CRISPRa kit - CRISPR gene activation of human paired box 3
Symbol:	PAX3
UniProt ID:	<a href="#">P23760</a>
Synonyms:	CDHS; HUP2; WS1; WS3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)