

Product datasheet for **GA103401**

Human PBX3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PBX3 CRISPRa kit - CRISPR gene activation of human PBX homeobox 3
Symbol:	PBX3
UniProt ID:	<u>P40426</u>
Synonyms:	OTTHUMP00000022145; OTTHUMP00000064213; pre-B-cell leukemia homeobox 3; pre-B-cell leukemia transcription factor 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)