

Product datasheet for **GA103404**

Human PCBP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCBP1 CRISPRa kit - CRISPR gene activation of human poly(rC) binding protein 1
Symbol:	PCBP1
UniProt ID:	Q15365
Synonyms:	HEL-S-85; hnRNP-E1; hnRNP-X; HNRPE1; HNRPX
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)