

Product datasheet for **GA103409**

Human PCDHGC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCDHGC3 CRISPRa kit - CRISPR gene activation of human protocadherin gamma subfamily C, 3
Symbol:	PCDHGC3
UniProt ID:	Q9UN70
Synonyms:	PC43; PCDH-GAMMA-C3; PCDH2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)