

Product datasheet for **GA103417**

Human PCMT1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCMT1 CRISPRa kit - CRISPR gene activation of human protein-L-isoaspartate (D-aspartate) O-methyltransferase
Symbol:	PCMT1
UniProt ID:	P22061
Synonyms:	PIMT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)