

Product datasheet for **GA103423**

Human PC1/3 (PCSK1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCSK1 CRISPRa kit - CRISPR gene activation of human proprotein convertase subtilisin/kexin type 1
Symbol:	PCSK1
UniProt ID:	P29120
Synonyms:	BMIQ12; NEC1; PC1; PC3; SPC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)