

Product datasheet for **GA103425**

Human Proprotein Convertase 2 (PCSK2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCSK2 CRISPRa kit - CRISPR gene activation of human proprotein convertase subtilisin/kexin type 2
Symbol:	PCSK2
UniProt ID:	P16519
Synonyms:	NEC-2; NEC 2; NEC2; PC2; SPC2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)