

## Product datasheet for **GA103431**

### Human PD1 (PDCD1) activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | PDCD1 CRISPRa kit - CRISPR gene activation of human programmed cell death 1 |
| Symbol:       | PDCD1   |
| UniProt ID:   | <a href="#">Q15116</a>  |
| Synonyms:     | CD279; hPD-1; hPD-I; hSLE1; PD-1; PD1; SLEB2                                |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)     |

#### Product images:



[View online >](#)