

## Product datasheet for **GA103437**

### Human PDE3B activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | PDE3B CRISPRa kit - CRISPR gene activation of human phosphodiesterase 3B |
| Symbol:       | PDE3B  |
| UniProt ID:   | <a href="#">Q13370</a>   |
| Synonyms:     | cGIPDE1; HcGIP1  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)  |

#### Product images:



[View online >](#)