


Product datasheet for **GA103467**

Human PDPK1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PDPK1 CRISPRa kit - CRISPR gene activation of human 3-phosphoinositide dependent protein kinase 1
Symbol:	PDPK1
UniProt ID:	<u>O15530</u>
Synonyms:	PDK1; PDPK2; PDPK2P; PRO0461
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)