

Product datasheet for **GA103468**

Human SLC26A4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC26A4 CRISPRa kit - CRISPR gene activation of human solute carrier family 26 member 4
Symbol:	SLC26A4
Components:	GA103468G1, SLC26A4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103468G2, SLC26A4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103468G3, SLC26A4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O43511</u>
Synonyms:	DFNB4; EVA; PDS; TDH2B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)