

Product datasheet for **GA103469**

Human ProDynorphin (PDYN) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PDYN CRISPRa kit - CRISPR gene activation of human prodynorphin
Symbol:	PDYN
UniProt ID:	P01213
Synonyms:	ADCA; PENKB; SCA23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)