

Product datasheet for **GA103471**

Human CD31 (PECAM1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PECAM1 CRISPRa kit - CRISPR gene activation of human platelet and endothelial cell adhesion molecule 1
Symbol:	PECAM1
UniProt ID:	P16284
Synonyms:	CD31; CD31/EndoCAM; endoCAM; GPIIA; PECA1; PECAM-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)