

Product datasheet for **GA103473**

Human PEG3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PEG3 CRISPRa kit - CRISPR gene activation of human paternally expressed 3
Symbol:	PEG3
UniProt ID:	<u>Q9GZU2</u>
Synonyms:	PW1; ZKSCAN22; ZNF904; ZSCAN24
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)