

Product datasheet for **GA103488**

Human PF4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PF4 CRISPRa kit - CRISPR gene activation of human platelet factor 4
Symbol:	PF4
UniProt ID:	<u>P02776</u>
Synonyms:	CXCL4; PF-4; SCYB4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)