

Product datasheet for **GA103494**

Human PFDN4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PFDN4 CRISPRa kit - CRISPR gene activation of human prefoldin subunit 4
Symbol:	PFDN4
Components:	GA103494G1, PFDN4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103494G2, PFDN4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103494G3, PFDN4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9NQ4
Synonyms:	C1; PFD4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)