

Product datasheet for **GA103495**

Human PFDN5 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PFDN5 CRISPRa kit - CRISPR gene activation of human prefoldin subunit 5 |
| Symbol: | PFDN5 |
| UniProt ID: | Q99471 |
| Synonyms: | MM-1; MM1; PFD5 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)