

Product datasheet for **GA103501**

Human PFKL activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PFKL CRISPRa kit - CRISPR gene activation of human phosphofructokinase, liver type
Symbol:	PFKL
UniProt ID:	P17858
Synonyms:	ATP-PFK; PFK-B; PFK-L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)