

Product datasheet for **GA103504**

Human PFKP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PFKP CRISPRa kit - CRISPR gene activation of human phosphofructokinase, platelet
Symbol:	PFKP
UniProt ID:	Q01813
Synonyms:	ATP-PFK; PFK-C; PFK-P; PFKF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

