

Product datasheet for **GA103511**

Human Pepsinogen II (PGC) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PGC CRISPRa kit - CRISPR gene activation of human progastricsin
Symbol:	PGC
UniProt ID:	<u>P20142</u>
Synonyms:	PEPC; PGII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)