

## Product datasheet for **GA103515**

### Human PGK1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PGK1 CRISPRa kit - CRISPR gene activation of human phosphoglycerate kinase 1
Symbol:	PGK1
UniProt ID:	<a href="#">P00558</a>
Synonyms:	HEL-S-68p; MIG10; PGKA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)