

Product datasheet for **GA103516**

Human PGK2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PGK2 CRISPRa kit - CRISPR gene activation of human phosphoglycerate kinase 2
Symbol:	PGK2
UniProt ID:	P07205
Synonyms:	dj417L20.2; HEL-S-272; PGKB; PGKPS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)