

## Product datasheet for **GA103517**

### Human PGM1 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | PGM1 CRISPRa kit - CRISPR gene activation of human phosphoglucomutase 1 |
| Symbol:       | PGM1  |
| UniProt ID:   | <a href="#">P36871</a>  |
| Synonyms:     | CDG1T; GSD14  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

#### Product images:



[View online >](#)