

Product datasheet for **GA103518**

Human PGM3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PGM3 CRISPRa kit - CRISPR gene activation of human phosphoglucosmutase 3
Symbol:	PGM3
UniProt ID:	O95394
Synonyms:	AGM1; IMD23; PAGM; PGM 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)