

## Product datasheet for **GA103530**

### Human Phosphorylase B (PHKB) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PHKB CRISPRa kit - CRISPR gene activation of human phosphorylase kinase regulatory subunit beta
Symbol:	PHKB
UniProt ID:	<a href="#">Q93100</a>
Synonyms:	DKFZp781E15103; FLJ41698
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)