

Product datasheet for **GA103535**

Human PI3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PI3 CRISPRa kit - CRISPR gene activation of human peptidase inhibitor 3
Symbol:	PI3
UniProt ID:	<u>P19957</u>
Synonyms:	cementoin; ESI; SKALP; WAP3; WFDC14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)