

## Product datasheet for **GA103538**

### Human SERPINB6 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | SERPINB6 CRISPRa kit - CRISPR gene activation of human serpin family B member 6 |
| Symbol:       | SERPINB6  |
| UniProt ID:   | <a href="#">P35237</a>  |
| Synonyms:     | CAP; DFNB91; MSTP057; PI-6; PI6; PTI; SPI3                                      |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)         |

#### Product images:



[View online >](#)