

Product datasheet for **GA103546**

Human PIGA activation kit by CRISPRa

Product data:

Product Type:

CRISPRa Kits

Description:

PIGA CRISPRa kit - CRISPR gene activation of human phosphatidylinositol glycan anchor biosynthesis class A

Symbol:

PIGA

UniProt ID:

[P37287](#)


Synonyms:

GPI3; MCAHS2; PIG-A; PNH1

Format:

3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)