

## Product datasheet for **GA103547**

### Human PIGC activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PIGC CRISPRa kit - CRISPR gene activation of human phosphatidylinositol glycan anchor biosynthesis class C
Symbol:	PIGC
UniProt ID:	<a href="#">Q92535</a>
Synonyms:	GPI2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)