

Product datasheet for **GA103550**

Human PIGR activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIGR CRISPRa kit - CRISPR gene activation of human polymeric immunoglobulin receptor
Symbol:	PIGR
UniProt ID:	P01833
Synonyms:	FLJ22667; MGC125361; MGC125362
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)