

Product datasheet for **GA103551**

Human PI 3 Kinase Class 2A (PIK3C2A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3C2A CRISPRa kit - CRISPR gene activation of human phosphatidylinositol-4-phosphate 3-kinase catalytic subunit type 2 alpha
Symbol:	PIK3C2A
UniProt ID:	<u>O00443</u>
Synonyms:	CPK; PI3-K-C2(ALPHA); PI3-K-C2A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)