

Product datasheet for **GA103552**

Human **PIK3C2B** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3C2B CRISPRa kit - CRISPR gene activation of human phosphatidylinositol-4-phosphate 3-kinase catalytic subunit type 2 beta
Symbol:	PIK3C2B
UniProt ID:	O00750
Synonyms:	C2-PI3K
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)