

## Product datasheet for **GA103553**

### Human **PIK3C2G** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	PIK3C2G CRISPRa kit - CRISPR gene activation of human phosphatidylinositol-4-phosphate 3-kinase catalytic subunit type 2 gamma
Symbol:	PIK3C2G
UniProt ID:	<u><a href="#">O75747</a></u>
Synonyms:	PI3K-C2-gamma; PI3K-C2GAMMA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)