

## Product datasheet for **GA103554**

### Human PI 3 Kinase Class 3 (PIK3C3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PIK3C3 CRISPRa kit - CRISPR gene activation of human phosphatidylinositol 3-kinase catalytic subunit type 3
Symbol:	PIK3C3
UniProt ID:	<a href="#">Q8NEB9</a>
Synonyms:	hVps34; VPS34; Vps34
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)