

## Product datasheet for **GA103555**

### Human PI 3 Kinase catalytic subunit alpha (PIK3CA) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PIK3CA CRISPRa kit - CRISPR gene activation of human phosphatidylinositol-4,5-bisphosphate 3-kinase catalytic subunit alpha
Symbol:	PIK3CA
UniProt ID:	<a href="#">P42336</a>
Synonyms:	CLOVE; CWS5; MCAP; MCM; MCMTc; p110-alpha; PI3K; PI3K-alpha
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)