

Product datasheet for **GA103556**

Human PIK3CB activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3CB CRISPRa kit - CRISPR gene activation of human phosphatidylinositol-4,5-bisphosphate 3-kinase catalytic subunit beta
Symbol:	PIK3CB
UniProt ID:	P42338
Synonyms:	P110BETA; PI3K; PI3KBETA; PIK3C1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)