

Product datasheet for **GA103559**

Human PI 3 Kinase catalytic subunit gamma (PIK3CG) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3CG CRISPRa kit - CRISPR gene activation of human phosphatidylinositol-4,5-bisphosphate 3-kinase catalytic subunit gamma
Symbol:	PIK3CG
UniProt ID:	P48736
Synonyms:	p110gamma; p120-PI3K; PI3CG; PI3K; PI3Kgamma; PIK3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)