

Product datasheet for **GA103561**

Human PI 3 Kinase p85 beta (PIK3R2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3R2 CRISPRa kit - CRISPR gene activation of human phosphoinositide-3-kinase regulatory subunit 2
Symbol:	PIK3R2
UniProt ID:	<u>O00459</u>
Synonyms:	MPPH; MPPH1; p85; p85-BETA; P85B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)