

Product datasheet for **GA103562**

Human PI4KA activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PI4KA CRISPRa kit - CRISPR gene activation of human phosphatidylinositol 4-kinase alpha
Symbol:	PI4KA
UniProt ID:	<u>P42356</u>
Synonyms:	PI4K-ALPHA; pi4K230; PIK4CA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)