

## Product datasheet for **GA103563**

### Human PI4KB activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PI4KB CRISPRa kit - CRISPR gene activation of human phosphatidylinositol 4-kinase beta
Symbol:	PI4KB
UniProt ID:	<u><a href="#">Q9UBF8</a></u>
Synonyms:	NPIK; PI4K-BETA; PI4K92; PI4KBETA; PI4KIIIBETA; PIK4CB
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)