

Product datasheet for **GA103565**

Human PIN4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIN4 CRISPRa kit - CRISPR gene activation of human peptidylprolyl cis/trans isomerase, NIMA-interacting 4
Symbol:	PIN4
UniProt ID:	Q9Y237
Synonyms:	EPVH; PAR14; PAR17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

