

Product datasheet for **GA103567**

Human PIP5K2 alpha (PIP4K2A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIP4K2A CRISPRa kit - CRISPR gene activation of human phosphatidylinositol-5-phosphate 4-kinase type 2 alpha
Symbol:	PIP4K2A
UniProt ID:	P48426
Synonyms:	PI5P4KA; PIP5K2A; PIP5KII-alpha; PIP5KIIA; PIPK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)