

Product datasheet for **GA103568**

Human PITPN (PITPNA) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PITPNA CRISPRa kit - CRISPR gene activation of human phosphatidylinositol transfer protein alpha
Symbol:	PITPNA
UniProt ID:	<u>Q00169</u>
Synonyms:	HEL-S-36; PI-TPalpha; PITPN; VIB1A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)