

Product datasheet for **GA103573**

Human Polycystin 2 (PKD2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PKD2 CRISPRa kit - CRISPR gene activation of human polycystin 2, transient receptor potential cation channel
Symbol:	PKD2
UniProt ID:	Q13563
Synonyms:	APKD2; Pc-2; PC2; PKD4; TRPP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)