

## Product datasheet for **GA103575**

### Human PKHD1 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | PKHD1 CRISPRa kit - CRISPR gene activation of human PKHD1 ciliary IPT domain containing fibrocystin/polyductin |
| Symbol:       | PKHD1  |
| UniProt ID:   | <a href="#">P08F94</a>   |
| Synonyms:     | ARPKD; FCYT; TIGM1   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)  |

#### Product images:



[View online >](#)